



Leopard Mask

Instructions:

1. Cut the mask out of orange paper. Cut the ear flaps on the solid lines and fold over along the dotted lines. Cut out the eyes.
2. Dip a sponge in black paint and make spot prints on the mask. (You may want to practice on a plain sheet of paper first.) Allow to dry.
3. Draw on the nose and mouth with a marker or crayon.
4. Cut the tongue out of pink paper. Tape it onto the back of the mask.
5. Thread white pipe cleaners through the holes shown on the pattern to make whiskers.
6. Fold the paper towel tube flat. Tape onto the back of the mask.

Tongue

